

Nitrite VACUettes® Kit

K-7004B/R-7002B: 0 - 300 ppm N

K-7004C/R-7002C: 0 - 3000 ppm N

Test Procedure

1. Fill the dilutor snapper cup to the -ml- mark with **distilled water** (fig. 1).
2. Fill the micro-test tube approximately halfway with the sample to be tested (fig. 2).
3. Place a VACUette tip firmly on to the ampoule tip.




Figure 1

4. Holding the VACUette almost horizontally, touch the tip to the contents of the microtest tube (fig. 2).




Figure 2

5. Place the VACUette between the vertical tip guides on the inside of the dilutor snapper cup. Snap the ampoule tip. The ampoule will fill leaving a bubble for mixing (fig. 3).




Figure 3

6. To mix the ampoule, invert it several times, allowing the bubble to travel from end to end.

7. Dry the ampoule. Obtain a test result **10 minutes** after snapping the tip.

8. Obtain a test result by placing the ampoule between the color standards until the best color match is found (fig 4).




Figure 4

NOTE: To convert to ppm nitrite (NO₂), multiply test result by 3.3.

Test Method

The Nitrite VACUettes®¹ test kit employs the azo dye formation method.^{2,3} In an acidic solution, nitrite diazotizes with a primary aromatic amine and then couples with another organic molecule to produce a highly colored azo dye. The resulting pink-orange color is proportional to the nitrite concentration in the sample.

1. VACUettes is a registered trademark of AquaPhoenix Scientific, LLC. U.S. Patent Nos. 4,537,747 & 4,596,780

2. APHA Standard Methods, 23rd ed., Method 4500-NO₂⁻ B - 2000

3. EPA Methods for Chemical Analysis of Water and Wastes, Method 354.1 (1983).

Safety Information

Read SDS before performing this test procedure. Wear safety glasses and protective gloves.

