

Multi-Analyte Photometers for Water Analysis

CHEMetrics' handheld, portable multi-analyte photometers are rugged and dependable solutions to your water testing needs in the field, the plant, or the laboratory. Intuitive and easy to use, these photometers allow for a quick and easy menu selection of 40+ pre-programmed analytes featuring the convenience of CHEMetrics® Vacu-vials® "snap and read" self-filling ampoules.

As new tests are available, a simple upload procedure from the CHEMetrics website updates the photometer with the most recent programs. Uploading takes only a few minutes and keeps the V-2000 and V-3000 current. The optional Power LabStation upgrades the portable V-3000 Photometer to a benchtop laboratory instrument.



Optional Accessories for V-3000

- A-0301** Data Management Software
- A-0302** Power LabStation
- A-0306** 28 mm cell with lid
- A-0307** RS232 to USB Adapter (for V-2000 & V-3000)

Specifications & Features

Features	V-2000 Photometer	V-3000 Photometer
Instrument Applicability	Portable	Portable / Benchtop
Display	LCD	Graphics / Backlit
Control Auto Shutoff	No	Yes
Power Supply Options	Battery	Battery Rechargeable Battery* Universal Cable / Plug*
Wavelengths (nm)	420, 520, 580, 610	436, 517, 557, 594, 610, 690
Data Interface Software	No	Yes
Cell Size	13 mm, 16 mm	13 mm, 28 mm
Language Selection	No	Yes: English, German, French, Spanish
Web-based Methods Update	Yes	Yes
Waterproof	IP67	IP67
Operating Temperature	0 to 45° C	0 to 50° C
Data Logging	100 points	100 points
Warranty	2 years	2 years

*Requires purchase of Power LabStation (A-0302)

CHEMetrics New V-0002 Verification Kit to check the performance of CHEMetrics® V-2000 Photometer.

See page 12 for more information.



Catalog No. Product Description

V-2000	Multi-Analyte Photometer
V-0002	V-2000 Verification Kit New!
V-3000	Multi-Analyte Photometer