

Updates to CHEMetrics' Product SDSs and Labels April 2024

Please be aware that CHEMetrics product SDSs and labels are actively being updated. This process began in December 2023 and will continue to be adjusted and implemented over the coming months. The goal of these changes is to ensure that all product regulatory information is compliant and clearly presented on our products. This means that there will be some new labels and in new locations.

Where applicable, labels with reagent information will be applied directly to the **BACK** of the individual trays that hold the ampoules (identified by an internal "WK-" or "WR-" product code). In addition to hazard information, the tray labels will include product lot identification and expiration date.




In products that include multiple reagent-containing components, such as Vacu-vials[®] and Titrets[®] Kits, each product component label will have full regulatory information. The kit label, applied to the **top panel of the kit box**, will carry a lot code identifier that captures lot number traceability for all components in the kit.

For CHEMetrics refills that do not contain multiple reagent-containing components, the label on the exterior of the product box will also have full regulatory lot and expiration date information.

Access to our SDSs will be changing as well. To obtain access to the new SDSs, either use the QR code on the product label, or go to the specific hyperlink on the product label (example: SDS: Available: sdsfetch.com/v1503ss). After using this path to get into SDS fetch, click on the "Test Kits" tab then on the Test kit ID to see SDSs for the components in your test kit. You can also then search for any CHEMetrics product number in the search bar.

Label Examples below (for K-1503 Ammonia Vacu-vials[®] Test Kit):

- Label on each tray of ampoules. The QR code directs you to the SDS for the ampoules.

 Ammonia Vacu-vials [®] Catalog No.: WK-1503	 <p>SIGNAL WORD: Danger HAZARD STATEMENT: Harmful if swallowed. May be harmful in contact with skin. Causes severe skin burns and eye damage. May cause an allergic skin reaction. May damage fertility. May damage the unborn child. May cause damage to organs through prolonged or repeated exposure. PRECAUTIONARY STATEMENT: Obtain, read and follow all safety instructions before use. Do not breathe dusts or mists. Wash hands thoroughly after handling. Do not touch eyes. Do not eat, drink or smoke when using this product. Contaminated work clothing should not be allowed out of the workplace. Wear eye protection/face protection. IF SWALLOWED: Get medical help. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN: Get medical help. IF ON SKIN: Wash with plenty of water. IF ON SKIN: Take off immediately all contaminated clothing. Immediately rinse with water for several minutes. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Immediately rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get emergency medical help immediately. IF exposed or concerned, get medical advice. Specific treatment (see on this label). If skin irritation or rash occurs: Get medical help. Take off contaminated clothing and wash it before reuse. Wash contaminated clothing before reuse. Store locked up. Dispose of contents/container to industrial combustion plant. For Laboratory Use Only. Not for Human Consumption.</p> <p>Supplied By: AquaPhoenix Scientific 860 Gitts Run Road, Hanover, PA 17331 Business Phone: (717) 632-1291 Emergency Phone: 800-255-3924 SDS Identifier: V1503SS SDS Available: http://sdsfetch.com/v1503ss</p> 
<p>Lot: 4GA0011 Expires: APR/24</p> <p>CONTAINS: Sodium Hydroxide CAS# 1310-73-2 Disodium Tetraborate Decahydrate CAS# 1303-96-4 Mercury II Iodide CAS# 7774-29-0 Potassium Iodide CAS# 7681-11-0</p>	

- Label on the bottle of Stabilizer Solution.


Stabilizer Solution
Catalog No.: S-1500

SIGNAL WORD: Warning
HAZARD STATEMENT: Causes mild skin irritation. Harmful to aquatic life.
PRECAUTIONARY STATEMENT: Avoid release to the environment. If skin irritation occurs: Get medical help. Dispose of contents/container to industrial combustion plant. For Laboratory Use Only. Not for Human Consumption.

CONTAINS:
Potassium Sodium Tartrate CAS# 6381-59-5
Potassium Hydroxide CAS# 1310-58-3
SDS Available: <http://sdsfetch.com/v1500ss>

Supplied By: AquaPhoenix Scientific 860 Gitts Run Road, Hanover, PA 17331 Lot: 4GA0012
Business #: (717) 632-1291 Emergency #: 800-255-3924 Expires: JUL/26
SDS Identifier: V1500SS

- Label on top panel of kit box (carries the kit catalog number for ordering purposes).



Supplied By: AquaPhoenix Scientific
860 Gitts Run Road, Hanover, PA 17331
Business Phone: (717) 632-1291
Emergency Phone: 800-255-3924

Lot: 4A036

If you have any questions or concerns about the labeling of our products, please contact our customer service online at <https://www.chemetrics.com/contact-us-customer-service/> or call us at 800.356.3072.